

Pressure Ulcers Made Simple



1 Assess

- Complete skin assessment: classify in accordance with EPUAP-NPUAP Pressure Ulcer Classification System¹.
- For categories 2-4 complete wound assessment:
 - Wound bed appearance (tissue type and %: slough, necrosis, granulation, suspected biofilm).
 - Size (length, width, depth, tunneling).
 - Exudate (colour, consistency, level).
 - Associated pain and/or odour.
 - Peri-wound skin condition (swelling, discolouration, maceration).
 - Signs and symptoms of infection (pain, odour, heat, redness, swelling, purulence).





2 Manage

- Cleanse, and where necessary debride, the wound to remove barriers to healing, e.g. slough, necrosis, biofilm.
- Select a dressing that¹:
 - Manages exudate
 - Protects peri-wound skin
 - Maintains a moist wound environment
 - Is comfortable for the patient
- Consider an anti-microbial dressing for ulcers that are infected or at risk of infection.

3 Monitor


- Re-assess the ulcer at each dressing change and document.
 - **Note:** Reverse grading should never be used to describe the healing of a pressure ulcer.²
- Review dressing selection considering changes in ulcer depth and symptoms.
- Consider the need to continue to protect skin integrity once the ulcer has healed.

PRESSURE ULCER CLASSIFICATION¹ & DRESSING SELECTION*


CATEGORY I	CATEGORY II	CATEGORY III	CATEGORY IV
			
<p>Intact skin with non-blanchable redness of a localized area usually over a bony prominence. Darkly pigmented skin may not have visible blanching.</p> <p>AQUACELTM Foam</p> <p>DuoDERM[®] Extra Thin</p>	<p>Partial thickness loss of dermis presenting as a shallow open ulcer with a red pink wound bed, without slough. May also present as an intact or open/ruptured serum-filled or sero-sanguinous filled blister.</p> <p>AQUACELTM Foam</p> <p>AQUACELTM Ag+ Foam[†]</p>	<p>Full thickness tissue loss. Subcutaneous fat may be visible but bone, tendon or muscle are not exposed. Slough may be present but does not obscure the depth of tissue loss. May include undermining and tunneling.</p> <p>AQUACELTM ExtraTM or AQUACELTM Ag+ ExtraTM† + AQUACELTM Foam</p>	<p>Full thickness tissue loss with exposed bone, tendon or muscle. Slough or eschar may be present. Often includes undermining and tunneling.</p> <p>AQUACELTM ExtraTM or AQUACELTM Ag+ ExtraTM† + AQUACELTM Foam</p> <p><small>AQUACELTM and AQUACELTM Ag+ Ribbon Dressings are available.</small></p>

The type of dressing may change over time as the ulcer heals or deteriorates.[‡]


DRESSING APPLICATION TIPS*



- AQUACELTM ExtraTM and AQUACELTM Ag+ ExtraTM dressings should overlap at least 1 cm onto the skin surrounding the wound.
- For cavity wounds AQUACELTM and AQUACELTM Ag+ Ribbon dressings are recommended.



- When dressing deep wounds only fill to 80% to allow for dressing expansion on contact with wound fluid.
- When using AQUACELTM and AQUACELTM Ag+ Ribbon dressings, leave 2.5 cm length of ribbon outside of the cavity to aid removal.



- Perfect Partners**
- Use AQUACELTM Foam as a cover dressing over AQUACELTM ExtraTM or AQUACELTM Ag+ Extra dressings.
 - The absorbent pad of AQUACELTM Foam dressings should overlap the wound by at least 1 cm.

1-800-XXX-XXXX / AQUACEL.com/foam

AQUACEL, AQUACEL Extra and DuoDERM are trademarks of ConvaTec Inc. and are registered trademarks in the US. ©2014 ConvaTec Inc. AP-014364-MM

* Please refer to pack insert for complete information on indications and usage for each product.

† Consider an anti-microbial dressing for ulcers that are infected or at risk of infection.

‡ Reverse grading should never be used to describe the healing of a pressure ulcer.

Reference: 1. European Pressure Ulcer Advisory Panel and National Pressure Ulcer Advisory Panel. Treatment of pressure ulcers: Quick Reference Guide. Washington DC: National Pressure Ulcer Advisory Panel; 2009. 2. http://www.npuap.org/wp-content/uploads/2012/01/NPUAP_position_on_staging-final.pdf

